

## GazeToMouse Toolchain

Generated by Doxygen 1.8.13



# Contents

<b>1</b>	<b>v2.3.0</b>	<b>1</b>
<b>2</b>	<b>Toolset to Control Tobii Eye Tracker 4C</b>	<b>5</b>
<b>3</b>	<b>Hierarchical Index</b>	<b>7</b>
3.1	Class Hierarchy . . . . .	7
<b>4</b>	<b>Class Index</b>	<b>9</b>
4.1	Class List . . . . .	9
<b>5</b>	<b>Class Documentation</b>	<b>11</b>
5.1	GazeHelper.ConfigItem Class Reference . . . . .	11
5.1.1	Detailed Description . . . . .	11
5.2	GazeHelper.EyeTrackerCore Class Reference . . . . .	12
5.2.1	Detailed Description . . . . .	13
5.2.2	Constructor & Destructor Documentation . . . . .	13
5.2.2.1	EyeTrackerCore() [1/2] . . . . .	13
5.2.2.2	EyeTrackerCore() [2/2] . . . . .	13
5.2.3	Member Function Documentation . . . . .	14
5.2.3.1	Dispose() . . . . .	14
5.3	GazeHelper.EyeTrackerPro Class Reference . . . . .	14
5.3.1	Detailed Description . . . . .	15
5.3.2	Constructor & Destructor Documentation . . . . .	15
5.3.2.1	EyeTrackerPro() . . . . .	15
5.3.3	Member Function Documentation . . . . .	16

5.3.3.1	IsLicenseOk()	16
5.4	GazeHelper.GazeDataArgs Class Reference	16
5.4.1	Detailed Description	17
5.4.2	Constructor & Destructor Documentation	17
5.4.2.1	GazeDataArgs() [1/2]	18
5.4.2.2	GazeDataArgs() [2/2]	18
5.5	GazeToMouse.GazeToMouse Class Reference	19
5.5.1	Detailed Description	20
5.5.2	Member Function Documentation	20
5.5.2.1	Main()	20
5.6	GazeHelper.JsonConfigParser Class Reference	20
5.6.1	Detailed Description	21
5.6.2	Constructor & Destructor Documentation	21
5.6.2.1	JsonConfigParser()	21
5.6.3	Member Function Documentation	21
5.6.3.1	GetDefaultConfig()	21
5.6.3.2	ParseJsonConfig()	22
5.6.3.3	SerializeJsonConfig()	22
5.7	GazeHelper.MouseHider Class Reference	22
5.7.1	Detailed Description	22
5.7.2	Constructor & Destructor Documentation	23
5.7.2.1	MouseHider()	23
5.7.3	Member Function Documentation	23
5.7.3.1	HideCursor()	23
5.7.3.2	ShowCursor()	23
5.8	GazeHelper.MouseTracker Class Reference	24
5.8.1	Detailed Description	25
5.8.2	Constructor & Destructor Documentation	25
5.8.2.1	MouseTracker()	25
5.8.3	Member Function Documentation	25

5.8.3.1	Dispose()	25
5.8.3.2	Start()	26
5.8.3.3	Stop()	26
5.9	GazeHelper.TrackerHandler Class Reference	26
5.9.1	Detailed Description	28
5.9.2	Constructor & Destructor Documentation	28
5.9.2.1	TrackerHandler()	28
5.9.3	Member Function Documentation	28
5.9.3.1	Dispose() [1/2]	28
5.9.3.2	Dispose() [2/2]	29
5.9.3.3	GazeDataHandler()	29
5.9.3.4	IsReady()	29
5.9.3.5	OnGazeDataReceived()	29
5.9.3.6	OnPropertyChanged()	30
5.9.3.7	OnTrackerDisabled()	30
5.9.3.8	OnTrackerDisabledTimeout()	30
5.9.3.9	OnTrackerEnabled()	31
5.9.4	Member Data Documentation	31
5.9.4.1	DeviceName	31
5.9.4.2	dialogBoxTimer	31
5.9.4.3	logger	31
5.9.4.4	trackerMessageBox	31
5.9.5	Property Documentation	32
5.9.5.1	State	32
5.9.6	Event Documentation	32
5.9.6.1	GazeDataReceived	32
5.9.6.2	PropertyChanged	32
5.9.6.3	TrackerDisabled	32
5.9.6.4	TrackerEnabled	32
5.10	GazeHelper.TrackerLogger Class Reference	33
5.10.1	Detailed Description	33
5.10.2	Constructor & Destructor Documentation	33
5.10.2.1	TrackerLogger()	33
5.10.3	Member Function Documentation	33
5.10.3.1	Debug()	33
5.10.3.2	DumpFatal()	34
5.10.3.3	Error()	34
5.10.3.4	Info()	34
5.10.3.5	Warning()	35
5.11	GazeHelper.TrackerMessageBox Class Reference	35
5.11.1	Detailed Description	36



# Chapter 1

## v2.3.0

### New Features

- A mouse tracker device can now be used instead of an eyetracker device. The mouse tracker logs the timestamp and the x and y coordinates of the mouse pointer whenever the mouse-move event is fired. The mouse tracker is used when the configuration field 'TrackerDevice' is set to the value 2.

### Improvements

- Rename the configuration field 'TobiiSDK' to 'TrackerDevice'.

## v2.2.0

### New Features

- Configuration file
  - Dump the configurations used for an experiment to a file at the "DataLogPath"
  - Allow to configure an experiment name which is used as a postfix of the dumped configuration file name
  - Consider the config file as invalid if not all required configuration parameters are defined
  - Consider the config file as invalid if unknown parameters are defined
  - Allow to configure whether to log data sets where all data is invalid (eyes closed, no subject in front of the screen, etc)
- Error Handling
  - Attach an error string to the output file, indicating errors that occurred during the run
  - Attach an error string to the dumped configuration file, indicating errors of the configuration

### Improvements

- Fall back to Core SDK if the license file cannot be applied to the device

## v2.1.0

### New Features

- Log eye origin coordinates
  - x, y, z coordinates of the left and the right eye
  - compute distance of the left and right eye to the eyetracker
  - compute the average distance of the two eyes

### Improvements

- Check the three format values and the column order individually to produce more specific log entries

## v2.0.1

### Bug Fix

- with SDK Pro, use system timestamp to cope with disconnected device
- fix the path in the z-tree sample file

## v2.0.0

### New Features

- Support for Tobii Pro SDK
  - apply license to eyetracker device at startup
  - logging of pupil diameter
  - logging of individual eye data
- Allow to configure column headers of output file

### Improvements

- Improved configuration options for the output file

## v1.0.0

### New Features

- Notify user with popup if eyetracker is not ready
- Allow to configure time interval for the software to wait for the eyetracker to become ready



#### Improvements

- Rename default output file for data from `<prefix>_data.txt` to `<prefix>_gaze.txt`

### v0.3.2

#### Improvements

- add header to the data log file.
- change the default value of allowed gaze data files.
- check and wait for ready state of the eye tracker before performing operations with it.

#### Bug Fix

- create a log file per machine to prevent concurrency conflicts.

### v0.3.1

#### Improvements

- ignore the option "HideMouse" when "ControlMouse" is disabled.

#### Bug Fix

- remove double log entry of mouse hiding and restoring event.

### v0.3.0

#### New Features

- allow to configure whether the gaze data is logged.
- allow to configure the maximum allowed amount of gaze data files in the output folder. Oldest files are deleted first.

#### Improvements

- limit the logfile size to 1MB. If the size is exceeded a new file is created. At any time only two log files are allowed, The older file is overwritten once both files exceed 1MB.

### v0.2.0

#### New Features

- allow to configure whether the mouse is controlled by the gaze of the subject or not.
- allow to configure the output format of the gaze data.

### v0.1.0

First release of the [GazeToMouse](#) toolset.

The toolset was tested on **Windows 7** in conjunction with **ztree v3.6.7** and Tobii **Eye Tracking Core v2.11.1.6952**.



## Chapter 2

# Toolset to Control Tobii Eye Tracker 4C

This repository contains the source code for multiple simple tools that allow to control a Tobii eye tracker from a 3rd party application. Specifically, this project aims at providing a set of executables that can be called from within [ztree](#) to allow eye tracker support for economic experiments. For more details please refer to the [documentation](#).

### Installation

The complete toolset package can be downloaded from the [release folder](#). The package contains the following executables:

- **TobiiCalibrate.exe** launches the Tobii calibration tool
- **TobiiGuestCalibrate.exe** launches the Tobii guest calibration
- **TobiiTest.exe** launches the Tobii test application
- **GazeToMouse.exe** enables the control of the mouse with the gaze of the user
- **GazeToMouseClose.exe** terminates `GazeToMouse.exe` gracefully
- **ShowMouse.exe** restores to mouse pointer if something went wrong

In order to run the executables the following files need to be placed in the same directory as the executables:

- GazeHelper.dll
- Newtonsoft.Json.dll
- Tobii.EyeX.Client.dll
- Tobii.Interaction.Model.dll
- Tobii.Interaction.Net.dll
- blank.cur
- config.json

Further, the Tobii engine must be running and the eye tracker must be enabled. To install the Tobii engine for the **Tobii Eye Tracker 4C** download the software from [here](#).

The Tobii engine starts the following service

```
Tobii.Service
```

and the following processes

```
Tobii EyeX Engine  
Tobii.EyeX.Interaction.exe  
Tobii.EyeX.Tray.exe
```

For quick starters, a simple **ztree sample program** is available.

## Release Notes

Information about the releases can be found in the [http://tpf.fluido.as:10012/TBI/TBI-tobii\\_eye\\_tracker\\_↔gaze/blob/master/CHANGELOG.md](http://tpf.fluido.as:10012/TBI/TBI-tobii_eye_tracker_↔gaze/blob/master/CHANGELOG.md) "CHANGELOG"

## Chapter 3

# Hierarchical Index

### 3.1 Class Hierarchy

This inheritance list is sorted roughly, but not completely, alphabetically:

Application	
GazeToMouse.GazeToMouse . . . . .	19
GazeHelper.ConfigItem . . . . .	11
EventArgs	
GazeHelper.GazeDataArgs . . . . .	16
IDisposable	
GazeHelper.TrackerHandler . . . . .	26
GazeHelper.EyeTrackerCore . . . . .	12
GazeHelper.EyeTrackerPro . . . . .	14
GazeHelper.MouseTracker . . . . .	24
INotifyPropertyChanged	
GazeHelper.TrackerHandler . . . . .	26
GazeHelper.JsonConfigParser . . . . .	20
GazeHelper.MouseHider . . . . .	22
GazeHelper.TrackerLogger . . . . .	33
Window	
GazeHelper.TrackerMessageBox . . . . .	35



## Chapter 4

# Class Index

### 4.1 Class List

Here are the classes, structs, unions and interfaces with brief descriptions:

<a href="#">GazeHelper.ConfigItem</a>	
configuration file class . . . . .	11
<a href="#">GazeHelper.EyeTrackerCore</a>	
Interface to the Tobii SDK Core engine . . . . .	12
<a href="#">GazeHelper.EyeTrackerPro</a>	
Interface to the Tobii SDK Pro engine . . . . .	14
<a href="#">GazeHelper.GazeDataArgs</a>	
The event argument class for Tobii eyetracker data . . . . .	16
<a href="#">GazeToMouse.GazeToMouse</a>	
Converts gaze data to mouse coordinates . . . . .	19
<a href="#">GazeHelper.JsonConfigParser</a>	
The config file "config.json" is parsed and its values are attributed to the <a href="#">ConfigItem</a> class. . . . .	20
<a href="#">GazeHelper.MouseHider</a>	
hide standard mouse pointer and restore it . . . . .	22
<a href="#">GazeHelper.MouseTracker</a>	
This class is used to hook into the system mouse events and track the position . . . . .	24
<a href="#">GazeHelper.TrackerHandler</a>	
The common interface for the Tobii eyetracker Engines Core and Pro . . . . .	26
<a href="#">GazeHelper.TrackerLogger</a>	
Simple logger class. . . . .	33
<a href="#">GazeHelper.TrackerMessageBox</a>	
Interaction logic for TrackerMessageBox.xaml . . . . .	35





## Chapter 5

# Class Documentation

### 5.1 GazeHelper.ConfigItem Class Reference

configuration file class

#### Properties

- string **ConfigName** [get, set]
- string **DataLogColumnOrder** [get, set]
- string[] **DataLogColumnTitle** [get, set]
- int **DataLogCount** [get, set]
- string **DataLogFormatDiameter** [get, set]
- string **DataLogFormatOrigin** [get, set]
- string **DataLogFormatTimeStamp** [get, set]
- bool **DataLogIgnoreInvalid** [get, set]
- string **DataLogPath** [get, set]
- bool **DataLogWriteOutput** [get, set]
- int **GazeFilterCore** [get, set]
- string **LicensePath** [get, set]
- bool **MouseControl** [get, set]
- bool **MouseHide** [get, set]
- string **MouseStandardIconPath** [get, set]
- int **ReadyTimer** [get, set]
- int **TrackerDevice** [get, set]
- string **TobiiApplicationPath** [get, set]
- string **TobiiCalibrate** [get, set]
- string **TobiiCalibrateArguments** [get, set]
- string **TobiiGuestCalibrate** [get, set]
- string **TobiiGuestCalibrateArguments** [get, set]
- string **TobiiTest** [get, set]

#### 5.1.1 Detailed Description

configuration file class


The documentation for this class was generated from the following file:

- GazeHelper/JsonConfigParser.cs


## 5.2 GazeHelper.EyeTrackerCore Class Reference

Interface to the Tobii SDK Core engine

Inheritance diagram for GazeHelper.EyeTrackerCore:



Collaboration diagram for GazeHelper.EyeTrackerCore:



### Public Member Functions

- [EyeTrackerCore](#) ([TrackerLogger](#) logger, int ready\_timer)  
*Initializes a new instance of the [EyeTrackerCore](#) class.*
- [EyeTrackerCore](#) ([TrackerLogger](#) logger, int ready\_timer, int gaze\_filter\_mode)  
*Initializes a new instance of the [EyeTrackerCore](#) class.*

### Protected Member Functions

- override void [Dispose](#) (bool disposing)  
*Releases unmanaged and - optionally - managed resources.*

## Additional Inherited Members

### 5.2.1 Detailed Description

Interface to the Tobii SDK Core engine

See also

[GazeHelper.TrackerHandler](#)

### 5.2.2 Constructor & Destructor Documentation

#### 5.2.2.1 EyeTrackerCore() [1/2]

```
GazeHelper.EyeTrackerCore.EyeTrackerCore (
    TrackerLogger logger,
    int ready_timer ) [inline]
```

Initializes a new instance of the [EyeTrackerCore](#) class.

Parameters

<i>logger</i>	The logger.
<i>ready_timer</i>	The ready timer.

#### 5.2.2.2 EyeTrackerCore() [2/2]

```
GazeHelper.EyeTrackerCore.EyeTrackerCore (
    TrackerLogger logger,
    int ready_timer,
    int gaze_filter_mode ) [inline]
```

Initializes a new instance of the [EyeTrackerCore](#) class.

Parameters

<i>logger</i>	The logger.
<i>ready_timer</i>	The ready timer.
<i>gaze_filter_mode</i>	The gaze filter mode.

## 5.2.3 Member Function Documentation

### 5.2.3.1 Dispose()

```
override void GazeHelper.EyeTrackerCore.Dispose (
    bool disposing ) [inline], [protected], [virtual]
```

Releases unmanaged and - optionally - managed resources.

#### Parameters

<i>disposing</i>	<code>true</code> to release both managed and unmanaged resources; <code>false</code> to release only unmanaged resources.
------------------	--

Reimplemented from [GazeHelper.TrackerHandler](#).


The documentation for this class was generated from the following file:

- GazeHelper/EyeTrackerCore.cs


## 5.3 GazeHelper.EyeTrackerPro Class Reference

Interface to the Tobii SDK Pro engine

Inheritance diagram for GazeHelper.EyeTrackerPro:



Collaboration diagram for GazeHelper.EyeTrackerPro:



## Public Member Functions

- [EyeTrackerPro](#) ([TrackerLogger](#) logger, int ready\_timer, string license\_path)  
*Initializes a new instance of the [EyeTrackerPro](#) class.*
- bool [IsLicenseOk](#) ()  
*Determines whether the license is applied to the eyetracker device*

## Additional Inherited Members

### 5.3.1 Detailed Description

Interface to the Tobii SDK Pro engine

See also

[GazeHelper.TrackerHandler](#)

### 5.3.2 Constructor & Destructor Documentation

#### 5.3.2.1 EyeTrackerPro()

```

GazeHelper.EyeTrackerPro.EyeTrackerPro (
    TrackerLogger logger,
    int ready_timer,
    string license_path ) [inline]
  
```

Initializes a new instance of the [EyeTrackerPro](#) class.

#### Parameters

<i>logger</i>	The logger.
<i>ready_timer</i>	The ready timer.
<i>license_path</i>	The license path.

### 5.3.3 Member Function Documentation

#### 5.3.3.1 IsLicenseOk()

```
bool GazeHelper.EyeTrackerPro.IsLicenseOk ( ) [inline]
```

Determines whether the license is applied to the eyetracker device

##### Returns

`true` if [is license ok]; otherwise, `false`.


The documentation for this class was generated from the following file:

- GazeHelper/EyeTrackerPro.cs


## 5.4 GazeHelper.GazeDataArgs Class Reference

The event argument class for Tobii eyetracker data

Inheritance diagram for GazeHelper.GazeDataArgs:



Collaboration diagram for GazeHelper.GazeDataArgs:



## Public Member Functions

- [GazeDataArgs](#) (TimeSpan timestamp, double xCoord, double yCoord)

*Initializes a new instance of the [GazeDataArgs](#) class.*

- [GazeDataArgs](#) (TimeSpan timestamp, double xCoord, double xCoordLeft, double xCoordRight, double yCoord, double yCoordLeft, double yCoordRight, bool isValidCoordLeft, bool isValidCoordRight, double dia, double diaLeft, double diaRight, bool isValidDiaLeft, bool isValidDiaRight, double xOriginLeft, double yOriginLeft, double zOriginLeft, double xOriginRight, double yOriginRight, double zOriginRight, double distOrigin, double distOriginLeft, double distOriginRight, bool isValidOriginLeft, bool isValidOriginRight)

*Initializes a new instance of the [GazeDataArgs](#) class.*

## Properties

- TimeSpan **Timestamp** [get]
- double **XCoord** [get]
- double **XCoordLeft** [get]
- double **XCoordRight** [get]
- double **YCoord** [get]
- double **YCoordLeft** [get]
- double **YCoordRight** [get]
- bool **IsValidCoordLeft** [get]
- bool **IsValidCoordRight** [get]
- double **Dia** [get]
- double **DiaLeft** [get]
- double **DiaRight** [get]
- bool **IsValidDiaLeft** [get]
- bool **IsValidDiaRight** [get]
- double **XOriginLeft** [get]
- double **XOriginRight** [get]
- double **YOriginLeft** [get]
- double **YOriginRight** [get]
- double **ZOriginLeft** [get]
- double **ZOriginRight** [get]
- double **DistOrigin** [get]
- double **DistOriginLeft** [get]
- double **DistOriginRight** [get]
- bool **IsValidOriginLeft** [get]
- bool **IsValidOriginRight** [get]

### 5.4.1 Detailed Description

The event argument class for Tobii eyetracker data

See also

System.EventArgs

### 5.4.2 Constructor & Destructor Documentation

### 5.4.2.1 GazeDataArgs() [1/2]

```
GazeHelper.GazeDataArgs.GazeDataArgs (
    TimeSpan timestamp,
    double xCoord,
    double yCoord ) [inline]
```

Initializes a new instance of the [GazeDataArgs](#) class.

#### Parameters

<i>timestamp</i>	The timestamp.
<i>xCoord</i>	The x coord of the gaze point.
<i>yCoord</i>	The y coord of the gaze point.

### 5.4.2.2 GazeDataArgs() [2/2]

```
GazeHelper.GazeDataArgs.GazeDataArgs (
    TimeSpan timestamp,
    double xCoord,
    double xCoordLeft,
    double xCoordRight,
    double yCoord,
    double yCoordLeft,
    double yCoordRight,
    bool isValidCoordLeft,
    bool isValidCoordRight,
    double dia,
    double diaLeft,
    double diaRight,
    bool isValidDiaLeft,
    bool isValidDiaRight,
    double xOriginLeft,
    double yOriginLeft,
    double zOriginLeft,
    double xOriginRight,
    double yOriginRight,
    double zOriginRight,
    double distOrigin,
    double distOriginLeft,
    double distOriginRight,
    bool isValidOriginLeft,
    bool isValidOriginRight ) [inline]
```

Initializes a new instance of the [GazeDataArgs](#) class.

#### Parameters

<i>timestamp</i>	The timestamp.
<i>xCoord</i>	The x coord of the gaze point.
<i>xCoordLeft</i>	The x coord of the gaze point of the left eye.
<i>xCoordRight</i>	The x coord of the gaze point of the right eye.



## Parameters

<i>yCoord</i>	The y coord of the gaze point.
<i>yCoordLeft</i>	The y coord of the gaze point of the left eye.
<i>yCoordRight</i>	The y coord of the gaze point of the right eye.
<i>isValidCoordLeft</i>	if set to <code>true</code> the gaze point coordinate of the left eye is valid.
<i>isValidCoordRight</i>	if set to <code>true</code> the gaze point coordinate of the right eye is valid.
<i>dia</i>	The average diameter of the pupils.
<i>diaLeft</i>	The diameter of the left pupil.
<i>diaRight</i>	The diameter of the right pupil.
<i>isValidDiaLeft</i>	if set to <code>true</code> the diameter of the left pupil is valid.
<i>isValidDiaRight</i>	if set to <code>true</code> the diameter of the right pupil is valid.
<i>xOriginLeft</i>	The x coord of the origin position of the left eye.
<i>yOriginLeft</i>	The y coord of the origin position of the left eye.
<i>zOriginLeft</i>	The z coord of the origin position of the left eye.
<i>xOriginRight</i>	The x coord of the origin position of the right eye.
<i>yOriginRight</i>	The y coord of the origin position of the right eye.
<i>zOriginRight</i>	The z coord of the origin position of the right eye.
<i>distOrigin</i>	The distance of the eye origin to the tracker.
<i>distOriginLeft</i>	The distance of the left eye origin to the tracker.
<i>distOriginRight</i>	The distance of the right eye origin to the tracker.
<i>isValidOriginLeft</i>	if set to <code>true</code> the origin point of the left eye is valid.
<i>isValidOriginRight</i>	if set to <code>true</code> the origin point of the right eye is valid.


The documentation for this class was generated from the following file:

- GazeHelper/TrackerHandler.cs


## 5.5 GazeToMouse.GazeToMouse Class Reference

Converts gaze data to mouse coordinates

Inheritance diagram for GazeToMouse.GazeToMouse:



Collaboration diagram for GazeToMouse.GazeToMouse:



### Static Public Member Functions

- static void [Main](#) ()  
*Defines the entry point of the application.*

#### 5.5.1 Detailed Description

Converts gaze data to mouse coordinates

See also

[System.Windows.Application](#)

#### 5.5.2 Member Function Documentation

##### 5.5.2.1 Main()

```
static void GazeToMouse.GazeToMouse.Main ( ) [inline], [static]
```

Defines the entry point of the application.

The documentation for this class was generated from the following file:

- [GazeToMouse/GazeToMouse.cs](#)

## 5.6 GazeHelper.JsonConfigParser Class Reference

The config file "config.json" is parsed and its values are attributed to the [ConfigItem](#) class.

## Public Member Functions

- [JsonConfigParser](#) ([TrackerLogger](#) logger)  
*Initializes a new instance of the [JsonConfigParser](#) class.*
- [ConfigItem ParseJsonConfig](#) ()  
*Parses the json configuration.*
- void [SerializeJsonConfig](#) ([ConfigItem](#) item, string path)  
*Serializes the json configuration object to a string and writes it to a file.*
- [ConfigItem GetDefaultConfig](#) ()  
*Gets the default configuration values.*

### 5.6.1 Detailed Description

The config file "config.json" is parsed and its values are attributed to the [ConfigItem](#) class.

### 5.6.2 Constructor & Destructor Documentation

#### 5.6.2.1 JsonConfigParser()

```
GazeHelper.JsonConfigParser.JsonConfigParser (  
    TrackerLogger logger ) [inline]
```

Initializes a new instance of the [JsonConfigParser](#) class.

#### Parameters

<i>logger</i>	The logger.
---------------	-------------

### 5.6.3 Member Function Documentation

#### 5.6.3.1 GetDefaultConfig()

```
ConfigItem GazeHelper.JsonConfigParser.GetDefaultConfig ( ) [inline]
```

Gets the default configuration values.

#### Returns

the default configuration values.

### 5.6.3.2 ParseJsonConfig()

`ConfigItem` GazeHelper.JsonConfigParser.ParseJsonConfig ( ) [inline]

Parses the json configuration.

#### Returns

the updated `ConfigItem` class.

### 5.6.3.3 SerializeJsonConfig()

```
void GazeHelper.JsonConfigParser.SerializeJsonConfig (
    ConfigItem item,
    string path ) [inline]
```

Serializes the json configuration object to a string and writes it to a file.

#### Parameters

<i>item</i>	The json configuration item.
<i>path</i>	The path where the file will be written.

The documentation for this class was generated from the following file:

- GazeHelper/JsonConfigParser.cs

## 5.7 GazeHelper.MouseHider Class Reference

hide standard mouse pointer and resore it

### Public Member Functions

- `MouseHider` (`TrackerLogger` logger)  
*Initializes a new instance of the `MouseHider` class.*
- void `HideCursor` ()  
*Hides the cursor.*
- void `ShowCursor` (string pathToCur)  
*Shows the cursor.*

### 5.7.1 Detailed Description

hide standard mouse pointer and resore it

## 5.7.2 Constructor & Destructor Documentation

### 5.7.2.1 MouseHider()

```
GazeHelper.MouseHider.MouseHider (
    TrackerLogger logger ) [inline]
```

Initializes a new instance of the [MouseHider](#) class.

#### Parameters

<i>logger</i>	The logger.
---------------	-------------

## 5.7.3 Member Function Documentation

### 5.7.3.1 HideCursor()

```
void GazeHelper.MouseHider.HideCursor ( ) [inline]
```

Hides the cursor.

Hides the standard mouse pointer by replacing the current icon with a transparent icon.

### 5.7.3.2 ShowCursor()

```
void GazeHelper.MouseHider.ShowCursor (
    string pathToCur ) [inline]
```

Shows the cursor.

the standard mouse pointer by replacing the current icon with the standard mouse pointer icon

#### Parameters

<i>pathToCur</i>	The path to the standard mouse pointer icon.
------------------	--


The documentation for this class was generated from the following file:

- GazeHelper/MouseHider.cs


## 5.8 GazeHelper.MouseTracker Class Reference

This class is used to hook into the system mouse events and track the position

Inheritance diagram for GazeHelper.MouseTracker:



Collaboration diagram for GazeHelper.MouseTracker:



### Public Member Functions

- [MouseTracker](#) ([TrackerLogger](#) logger, int ready\_timer)  
*Initializes a new instance of the [MouseTracker](#) class.*
- void [Start](#) ()  
*Hooks the callback function `HookCallback(int, IntPtr, IntPtr)` to mouse events.*
- void [Stop](#) ()  
*Removes to mouse event hook.*

### Protected Member Functions

- override void [Dispose](#) (bool disposing)  
*Releases unmanaged and - optionally - managed resources.*

## Additional Inherited Members

### 5.8.1 Detailed Description

This class is used to hook into the system mouse events and track the position

See also

[GazeHelper.TrackerHandler](#)

### 5.8.2 Constructor & Destructor Documentation

#### 5.8.2.1 MouseTracker()

```
GazeHelper.MouseTracker.MouseTracker (
    TrackerLogger logger,
    int ready_timer ) [inline]
```

Initializes a new instance of the [MouseTracker](#) class.

Parameters

<i>logger</i>	The logger.
<i>ready_timer</i>	The ready timer.

### 5.8.3 Member Function Documentation

#### 5.8.3.1 Dispose()

```
override void GazeHelper.MouseTracker.Dispose (
    bool disposing ) [inline], [protected], [virtual]
```

Releases unmanaged and - optionally - managed resources.

Parameters

<i>disposing</i>	true to release both managed and unmanaged resources; false to release only unmanaged resources.
------------------	--

Reimplemented from [GazeHelper.TrackerHandler](#).

### 5.8.3.2 Start()

```
void GazeHelper.MouseTracker.Start ( ) [inline]
```

Hooks the callback function HookCallback(int, IntPtr, IntPtr) to mouse events.

### 5.8.3.3 Stop()

```
void GazeHelper.MouseTracker.Stop ( ) [inline]
```

Removes to mouse event hook.


The documentation for this class was generated from the following file:

- GazeHelper/MouseTracker.cs


## 5.9 GazeHelper.TrackerHandler Class Reference

The common interface for the Tobii eyetracker Engines Core and Pro

Inheritance diagram for GazeHelper.TrackerHandler:



Collaboration diagram for GazeHelper.TrackerHandler:





## Public Member Functions

- delegate void [GazeDataHandler](#) (Object sender, [GazeDataArgs](#) e)  
*Event handler for gaze data events of the eyetracker*
- [TrackerHandler](#) ([TrackerLogger](#) logger, int ready\_timer, string device\_name)  
*Initializes a new instance of the EyeTrackerHandler class.*
- void [Dispose](#) ()  
*Performs application-defined tasks associated with freeing, releasing, or resetting unmanaged resources.*

## Public Attributes

- readonly string [DeviceName](#)  
*The name of the tracker device*

## Protected Member Functions

- virtual void [Dispose](#) (bool disposing)  
*Releases unmanaged and - optionally - managed resources.*
- bool [IsReady](#) ()  
*Determines whether this eye tracker is ready.*
- virtual void [OnGazeDataReceived](#) ([GazeDataArgs](#) e)  
*Called when [gaze data received].*
- virtual void [OnPropertyChanged](#) (string property\_name)  
*Called when when the state property of EyeTracker is changing.*
- virtual void [OnTrackerDisabled](#) (EventArgs e)  
*Raises the E:TrackerDisabled event.*
- void [OnTrackerDisabledTimeout](#) (Object source, System.Timers.ElapsedEventArgs e)  
*Called after a specified amount of time of the eyetracker not being ready.*
- virtual void [OnTrackerEnabled](#) (EventArgs e)  
*Raises the E:TrackerEnabled event.*

## Protected Attributes

- Timer [dialogBoxTimer](#)  
*Timer to control the apperance of the dialog box*
- [TrackerLogger](#) logger  
*The logger*
- [TrackerMessageBox](#) trackerMessageBox  
*The dialog box taht is controlled by the dialogBoxTimer*

## Properties

- EyeTrackingDeviceStatus [State](#) [get, set]  
*Gets or sets the state of the eyetracker device.*

## Events

- EventHandler [TrackerEnabled](#)  
*Occurs when [tracker enabled].*
- EventHandler [TrackerDisabled](#)  
*Occurs when [tracker disabled].*
- PropertyChangedEventHandler [PropertyChanged](#)  
*Occurs when a property value changes.*
- [GazeDataHandler](#) [GazeDataReceived](#)  
*Occurs when [gaze data received].*

### 5.9.1 Detailed Description

The common interface for the Tobii eyetracker Engines Core and Pro

#### See also

[System.ComponentModel.INotifyPropertyChanged](#), [System.IDisposable](#)

### 5.9.2 Constructor & Destructor Documentation

#### 5.9.2.1 TrackerHandler()

```
GazeHelper.TrackerHandler.TrackerHandler (
    TrackerLogger logger,
    int ready_timer,
    string device_name ) [inline]
```

Initializes a new instance of the EyeTrackerHandler class.

#### Parameters

<i>logger</i>	The logger.
<i>ready_timer</i>	The ready timer.
<i>device_name</i>	Name of the device.

### 5.9.3 Member Function Documentation

#### 5.9.3.1 Dispose() [1/2]

```
void GazeHelper.TrackerHandler.Dispose ( ) [inline]
```

Performs application-defined tasks associated with freeing, releasing, or resetting unmanaged resources.

**5.9.3.2 Dispose()** [2/2]

```
virtual void GazeHelper.TrackerHandler.Dispose (
    bool disposing ) [inline], [protected], [virtual]
```

Releases unmanaged and - optionally - managed resources.

**Parameters**

<i>disposing</i>	true to release both managed and unmanaged resources; false to release only unmanaged resources.
------------------	--

Reimplemented in [GazeHelper.MouseTracker](#), and [GazeHelper.EyeTrackerCore](#).

**5.9.3.3 GazeDataHandler()**

```
delegate void GazeHelper.TrackerHandler.GazeDataHandler (
    Object sender,
    GazeDataArgs e )
```

Event handler for gaze data events of the eyetracker

**Parameters**

<i>sender</i>	The sender.
<i>e</i>	The e.

**5.9.3.4 IsReady()**

```
bool GazeHelper.TrackerHandler.IsReady ( ) [inline], [protected]
```

Determines whether this eye tracker is ready.

**Returns**

true if this instance is ready; otherwise, false.

**5.9.3.5 OnGazeDataReceived()**

```
virtual void GazeHelper.TrackerHandler.OnGazeDataReceived (
    GazeDataArgs e ) [inline], [protected], [virtual]
```

Called when [gaze data received].

**Parameters**

<i>e</i>	The gaze data event data.
----------	---------------------------

**5.9.3.6 OnPropertyChanged()**

```
virtual void GazeHelper.TrackerHandler.OnPropertyChanged (
    string property_name ) [inline], [protected], [virtual]
```

Called when when the state property of EyeTracker is changing.

**Parameters**

<i>property_name</i>	Name of the property in WPF.
----------------------	------------------------------

**5.9.3.7 OnTrackerDisabled()**

```
virtual void GazeHelper.TrackerHandler.OnTrackerDisabled (
    EventArgs e ) [inline], [protected], [virtual]
```

Raises the E:TrackerDisabled event.

**Parameters**

<i>e</i>	The EventArgs instance containing the event data.
----------	---

**5.9.3.8 OnTrackerDisabledTimeout()**

```
void GazeHelper.TrackerHandler.OnTrackerDisabledTimeout (
    Object source,
    System.Timers.ElapsedEventArgs e ) [inline], [protected]
```

Called after a specified amount of time of the eyetracker not being ready.

**Parameters**

<i>source</i>	The source.
<i>e</i>	The System.Timers.ElapsedEventArgs instance containing the event data.

#### 5.9.3.9 OnTrackerEnabled()

```
virtual void GazeHelper.TrackerHandler.OnTrackerEnabled (
    EventArgs e ) [inline], [protected], [virtual]
```

Raises the E:TrackerEnabled event.

##### Parameters

<i>e</i>	The EventArgs instance containing the event data.
----------	---

### 5.9.4 Member Data Documentation

#### 5.9.4.1 DeviceName

```
readonly string GazeHelper.TrackerHandler.DeviceName
```

The name of the tracker device

#### 5.9.4.2 dialogBoxTimer

```
Timer GazeHelper.TrackerHandler.dialogBoxTimer [protected]
```

Timer to control the apperance of the dialog box

#### 5.9.4.3 logger

```
TrackerLogger GazeHelper.TrackerHandler.logger [protected]
```

The logger

#### 5.9.4.4 trackerMessageBox

```
TrackerMessageBox GazeHelper.TrackerHandler.trackerMessageBox [protected]
```

The dialog box taht is controlled by the dialogBoxTimer

## 5.9.5 Property Documentation

### 5.9.5.1 State

`EyeTrackingDeviceStatus GazeHelper.TrackerHandler.State [get], [set]`

Gets or sets the state of the eyetracker device.

The state.

## 5.9.6 Event Documentation

### 5.9.6.1 GazeDataReceived

`GazeDataHandler GazeHelper.TrackerHandler.GazeDataReceived`

Occurs when [gaze data received].

### 5.9.6.2 PropertyChanged

`PropertyChangedEventHandler GazeHelper.TrackerHandler.PropertyChanged`

Occurs when a property value changes.

### 5.9.6.3 TrackerDisabled

`EventHandler GazeHelper.TrackerHandler.TrackerDisabled`

Occurs when [tracker disabled].

### 5.9.6.4 TrackerEnabled

`EventHandler GazeHelper.TrackerHandler.TrackerEnabled`

Occurs when [tracker enabled].

The documentation for this class was generated from the following file:

- GazeHelper/TrackerHandler.cs

## 5.10 GazeHelper.TrackerLogger Class Reference

Simple logger class.

### Public Member Functions

- [TrackerLogger](#) ()  
*Initializes a new instance of the [TrackerLogger](#) class.*
- void [DumpFatal](#) (Exception e)  
*Dumps exception to a new file if it is not possible to write to the main log file.*
- void [Debug](#) (string message)  
*wrapper function for debug level logging.*
- void [Info](#) (string message)  
*wrapper function for info level logging*
- void [Warning](#) (string message)  
*wrapper function for warning level logging*
- void [Error](#) (string message)  
*wrapper function for error level logging*

### 5.10.1 Detailed Description

Simple logger class.

### 5.10.2 Constructor & Destructor Documentation

#### 5.10.2.1 TrackerLogger()

```
GazeHelper.TrackerLogger.TrackerLogger ( ) [inline]
```

Initializes a new instance of the [TrackerLogger](#) class.

### 5.10.3 Member Function Documentation

#### 5.10.3.1 Debug()

```
void GazeHelper.TrackerLogger.Debug (
    string message ) [inline]
```

wrapper function for debug level logging.

**Parameters**

<i>message</i>	The message.
----------------	--------------

**5.10.3.2 DumpFatal()**

```
void GazeHelper.TrackerLogger.DumpFatal (
    Exception e ) [inline]
```

Dumps exception to a new file if it is not possible to write to the main log file.

**Parameters**

<i>e</i>	The exception.
----------	----------------

**5.10.3.3 Error()**

```
void GazeHelper.TrackerLogger.Error (
    string message ) [inline]
```

wrapper function for error level logging

**Parameters**

<i>message</i>	The message.
----------------	--------------

**5.10.3.4 Info()**

```
void GazeHelper.TrackerLogger.Info (
    string message ) [inline]
```

wrapper function for info level logging

**Parameters**

<i>message</i>	The message.
----------------	--------------



#### 5.10.3.5 Warning()

```
void GazeHelper.TrackerLogger.Warning (
    string message ) [inline]
```

wrapper function for warning level logging

##### Parameters

<i>message</i>	The message.
----------------	--------------


The documentation for this class was generated from the following file:

- GazeHelper/Logger.cs


## 5.11 GazeHelper.TrackerMessageBox Class Reference

Interaction logic for TrackerMessageBox.xaml

Inheritance diagram for GazeHelper.TrackerMessageBox:



Collaboration diagram for GazeHelper.TrackerMessageBox:



### 5.11.1 Detailed Description

Interaction logic for TrackerMessageBox.xaml

The documentation for this class was generated from the following file:

- GazeHelper/TrackerMessageBox.xaml.cs

# Index

- Debug
  - GazeHelper::TrackerLogger, [33](#)
- DeviceName
  - GazeHelper::TrackerHandler, [31](#)
- dialogBoxTimer
  - GazeHelper::TrackerHandler, [31](#)
- Dispose
  - GazeHelper::EyeTrackerCore, [14](#)
  - GazeHelper::MouseTracker, [25](#)
  - GazeHelper::TrackerHandler, [28](#)
- DumpFatal
  - GazeHelper::TrackerLogger, [34](#)
- Error
  - GazeHelper::TrackerLogger, [34](#)
- EyeTrackerCore
  - GazeHelper::EyeTrackerCore, [13](#)
- EyeTrackerPro
  - GazeHelper::EyeTrackerPro, [15](#)
- GazeDataArgs
  - GazeHelper::GazeDataArgs, [17](#), [18](#)
- GazeDataHandler
  - GazeHelper::TrackerHandler, [29](#)
- GazeDataReceived
  - GazeHelper::TrackerHandler, [32](#)
- GazeHelper.ConfigItem, [11](#)
- GazeHelper.EyeTrackerCore, [12](#)
- GazeHelper.EyeTrackerPro, [14](#)
- GazeHelper.GazeDataArgs, [16](#)
- GazeHelper.JsonConfigParser, [20](#)
- GazeHelper.MouseHider, [22](#)
- GazeHelper.MouseTracker, [24](#)
- GazeHelper.TrackerHandler, [26](#)
- GazeHelper.TrackerLogger, [33](#)
- GazeHelper.TrackerMessageBox, [35](#)
- GazeHelper::EyeTrackerCore
  - Dispose, [14](#)
  - EyeTrackerCore, [13](#)
- GazeHelper::EyeTrackerPro
  - EyeTrackerPro, [15](#)
  - IsLicenseOk, [16](#)
- GazeHelper::GazeDataArgs
  - GazeDataArgs, [17](#), [18](#)
- GazeHelper::JsonConfigParser
  - GetDefaultConfig, [21](#)
  - JsonConfigParser, [21](#)
  - ParseJsonConfig, [21](#)
  - SerializeJsonConfig, [22](#)
- GazeHelper::MouseHider
  - HideCursor, [23](#)
  - MouseHider, [23](#)
  - ShowCursor, [23](#)
- GazeHelper::MouseTracker
  - Dispose, [25](#)
  - MouseTracker, [25](#)
  - Start, [25](#)
  - Stop, [26](#)
- GazeHelper::TrackerHandler
  - DeviceName, [31](#)
  - dialogBoxTimer, [31](#)
  - Dispose, [28](#)
  - GazeDataHandler, [29](#)
  - GazeDataReceived, [32](#)
  - IsReady, [29](#)
  - logger, [31](#)
  - OnGazeDataReceived, [29](#)
  - OnPropertyChanged, [30](#)
  - OnTrackerDisabled, [30](#)
  - OnTrackerDisabledTimeout, [30](#)
  - OnTrackerEnabled, [30](#)
  - PropertyChanged, [32](#)
  - State, [32](#)
  - TrackerDisabled, [32](#)
  - TrackerEnabled, [32](#)
  - TrackerHandler, [28](#)
  - trackerMessageBox, [31](#)
- GazeHelper::TrackerLogger
  - Debug, [33](#)
  - DumpFatal, [34](#)
  - Error, [34](#)
  - Info, [34](#)
  - TrackerLogger, [33](#)
  - Warning, [34](#)
- GazeToMouse.GazeToMouse, [19](#)
- GazeToMouse::GazeToMouse
  - Main, [20](#)
- GetDefaultConfig
  - GazeHelper::JsonConfigParser, [21](#)
- HideCursor
  - GazeHelper::MouseHider, [23](#)
- Info
  - GazeHelper::TrackerLogger, [34](#)
- IsLicenseOk
  - GazeHelper::EyeTrackerPro, [16](#)
- IsReady
  - GazeHelper::TrackerHandler, [29](#)

- JsonConfigParser
  - GazeHelper::JsonConfigParser, [21](#)
- logger
  - GazeHelper::TrackerHandler, [31](#)
- Main
  - GazeToMouse::GazeToMouse, [20](#)
- MouseHider
  - GazeHelper::MouseHider, [23](#)
- MouseTracker
  - GazeHelper::MouseTracker, [25](#)
- OnGazeDataReceived
  - GazeHelper::TrackerHandler, [29](#)
- OnPropertyChanged
  - GazeHelper::TrackerHandler, [30](#)
- OnTrackerDisabled
  - GazeHelper::TrackerHandler, [30](#)
- OnTrackerDisabledTimeout
  - GazeHelper::TrackerHandler, [30](#)
- OnTrackerEnabled
  - GazeHelper::TrackerHandler, [30](#)
- ParseJsonConfig
  - GazeHelper::JsonConfigParser, [21](#)
- PropertyChanged
  - GazeHelper::TrackerHandler, [32](#)
- SerializeJsonConfig
  - GazeHelper::JsonConfigParser, [22](#)
- ShowCursor
  - GazeHelper::MouseHider, [23](#)
- Start
  - GazeHelper::MouseTracker, [25](#)
- State
  - GazeHelper::TrackerHandler, [32](#)
- Stop
  - GazeHelper::MouseTracker, [26](#)
- TrackerDisabled
  - GazeHelper::TrackerHandler, [32](#)
- TrackerEnabled
  - GazeHelper::TrackerHandler, [32](#)
- TrackerHandler
  - GazeHelper::TrackerHandler, [28](#)
- TrackerLogger
  - GazeHelper::TrackerLogger, [33](#)
- trackerMessageBox
  - GazeHelper::TrackerHandler, [31](#)
- Warning
  - GazeHelper::TrackerLogger, [34](#)